




## Barkey S-line

### Safety

- Integrated temperature sensors
- Overtemperature alarm
- Smooth surfaces, therefore simple to wash and disinfect
- Operational safety thanks to preset target temperature value

### Prevention of hypothermia

Warmth is the key to patients' wellbeing. Barkey's inexpensive infusion and blood warmer S-line helps to ease the treatment for all patients and accelerates the heating process. As a result care personnel has less work and treatment costs are reduced.



Product	Barkey S-line
Mains supply	100-240 V, 50/60 Hz
Power consumption	115-139 VA
Set temperature	+37°C oder +39°C
Length warming sleeve	1.5 m
For I.V. set diameter	4.0–5.0 mm or 6.0–7.0 mm
Access line	3 m
Weight	approx. 1.3 kg

MDD classification	Ila
Protection class	I, protective earthing
Moisture protection	IP X2
Degree of protection	BF, defi-proof / internal and external treatment to the patient
Certification	CE 0123
Barkey is certified	DIN EN ISO 13485