

## *Introducing the NEW IMPROVED*

**MERCURY**  
MEDICAL®

### *Disposable Manometer*

*Now color-coded for improved visibility of pressure ranges.*



Color-Coded Manometer: green 0 - 20 cm H<sub>2</sub>O,  
yellow between 20 - 40 cm H<sub>2</sub>O and red over 40 cm H<sub>2</sub>O.



The Mercury Medical Disposable Pressure Manometer is the first disposable manometer to reliably monitor BOTH proper inflation pressure and PEEP pressure during manual ventilation with Mercury's CPR and Hyperinflation Systems. And now the Mercury Medical Disposable Pressure Manometer is available with a color-coded label for improved visibility of pressure ranges.

**MERCURY**  
MEDICAL®

### *Disposable Manometer...*

*Ready when you need it and Now easier to read!*



# New Improved Disposable Manometer

measures BOTH Airway and PEEP pressures and is now color-coded

The Mercury Medical Disposable Color-Coded Manometer allows for in-line monitoring of BOTH airway pressure and PEEP pressures, when used in conjunction with Mercury's CPR resuscitation bags.

The Mercury Medical Disposable Color-Coded Manometer is ready when you need a manometer so you don't need to search or deal with a broken reusable manometer again!

**Only the Mercury Medical Color-Coded Manometer offers all the following:**

- New color-coding for improved visibility of pressure ranges.
- In-line use attaches directly onto CPR and hyperinflation resuscitation bags reducing the need to look away from the patient while monitoring pressure.\*
- Monitoring of BOTH airway and PEEP pressures.
- Disposable feature to address potential infection control issues and eliminate search for misplaced or broken reusable devices.
- Latex-Free
- Accuracy of  $\pm 3$  cm H<sub>2</sub>O up to 15 cm H<sub>2</sub>O
- Accuracy of  $\pm 5$  cm H<sub>2</sub>O over 15 cm H<sub>2</sub>O



Squeeze the resuscitation bag to ventilate and the manometer verifies airway pressure.



Release resuscitation bag squeeze on exhalation and the manometer verifies PEEP pressure (when PEEP valve is attached to the resuscitator).

## Ordering Information

Item	Description	Quantity
10-55350	Disposable Pressure Manometer, 0-60 cm H <sub>2</sub> O	20/box
10-55351	Disposable Pressure Manometer with stand and 3' tubing, 0-60 cm H <sub>2</sub> O	12/box
10-55356	Manometer with Tee, Left Orientation, 0-60 cm H <sub>2</sub> O	20/box
10-55357	Manometer with Tee, Right Orientation, 0-60 cm H <sub>2</sub> O	20/box



**Disposable Manometer... Ready when you need it and Now easier to read!**

*\*Many Mercury CPR and Hyperinflation resuscitators come equipped with in-line, color-coded manometers.*

Distributed by:



P O Box 2541 Southport, QLD 4215  
Ph: 1300-786-454 F: 5591 8522  
E; sales@questsurgical.com.au  
ABN 33 118 090 618 ACN 118 090 618

Rev. A