

Manometer *(Mercury Medical Disposable Color-Coded Manometer)*

The Disposable Pressure Manometer is the first disposable manometer to reliably monitor BOTH proper inflation pressure and PEEP pressure during manual ventilation with Mercury's CPR and Hyperinflation Systems. And now, the Disposable Pressure Manometer is available with a color-coded label for improved visibility of pressure ranges.



Color-Coded Manometer: green 0-20 cm H₂O, yellow, between 20-40 cm H₂O and red over 40 cm H₂O

Features

- The Mercury Medical Disposable Color-Coded Manometer allows for in-line monitoring of both airway pressure and PEEP pressures, when used in conjunction with Mercury's CPR and Hyperinflation resuscitation bags.
- The Mercury Medical Disposable Color-Coded Manometer is ready when you need a manometer so you don't need to search or deal with a broken reusable manometer again!
- Disposable
- Latex-Free

Benefits

- New color-coding for improved visibility of pressure ranges.
- In-line use attaches directly onto CPR and Hyperinflation resuscitation bags reducing the need to look away from the patient while monitoring pressure.
- Monitoring of both airway and PEEP pressures.
- Accuracy of ± 3 cm H₂O up to 15 cm H₂O.
- Accuracy of ± 5 cm H₂O over 15 cm H₂O.

ORDERING INFORMATION

Part#	Description	Quantity
#10-55350	Disposable Pressure Manometer, 0-60 cm H ₂ O	20/Box
#10-55351	Disposable Pressure Manometer with stand and 3' tubing, 0-60 cm H ₂ O	12/Box
#10-55352	Disposable Pressure Manometer, 0-40 cm H ₂ O	20/Box
#10-55353	Manometer with Tee, Left Orientation, 0-40 cm H ₂ O	20/Box
#10-55354	Manometer with Tee, Right Orientation, 0-40 cm H ₂ O	20/Box
#10-55356	Manometer with Tee, Left Orientation, 0-60 cm H ₂ O	20/Box
#10-55357	Manometer with Tee, Right Orientation, 0-60 cm H ₂ O	20/Box