

SURGICAL GOWNS


CLASSIC AND COOLMAX RANGES

Surgical Gowns with varying degrees of protection.

ALL OF OUR CLASSIC & COOLMAX GOWNS:

- are T-cut with velcro neckline tabs for easy adjustment
- have high quality cuffs made with a soft cotton fabric
- meet AAMI according to Level 3 or 4 surgical gown standards

- come wrapped in a 91cm x 91cm SMS field with 40cm x 68cm huck towel
- are sterile and latex free
- come in cartons of 25.



AAMI LEVEL 3

AAMI LEVEL 4

	CLASSIC UNREINFORCED				CLASSIC POLY REINFORCED				COOLMAX SURGICAL GOWN			
SIZE	M	L	XL	XLXL	M	L	XL	XLXL	M	L	XL	XLXL
CODE	DEF5000A	DEF5001A	DEF5002A	DEF5004A	DEF5005A	DEF5006A	DEF5007A	DEF5009A	DEF5010A	DEF5011A	DEF5012A	DEF5013A
FEATURES	<ul style="list-style-type: none"> • Comfortable and light • Soft material and cool to wear 				<ul style="list-style-type: none"> • Poly reinforcement on front panel and sleeves • Durable and impervious to fluid on panels • Comfortable to wear 				<ul style="list-style-type: none"> • Trilaminare Breathable Viral Barrier (BVB) material • High-tech breathable fabric which is fully impervious • Comfortable and cool to wear 			
THICKNESS	50 gsm				50 gsm + PE Zoned				68 gsm			
PROCEDURE	Dry and short				Medium or moderate fluid				High-fluid or lengthy			

MATERIALS KEY

- Unreinforced
- Poly Reinforced
- Trilaminare BVB

The clear choice for personal protection

SPEYES™ EYE SHIELDZ™



Specifications

- Length: 23cm
- Width: 30cm
- Vented foam with strap
- Peel away protective cover with tabs to prevent scratches
- Anti-fog
- Non-Sterile, Single use
- Latex Free
- ARTG 334977
- Meets AS/ NZS 1337

Shield your eyes from blood, bodily fluids or other risky splashes with quality eye protection.

The DeRoyal SPEyes™ Eye Shieldz™ range are designed to fit snugly for complete splash protection while UltraClear™ technology prevents glare, fog and static. Clip on and off designs enable quick change and minimise handling and contamination risk.

Give your healthcare providers a clear view to focus on what they do best.

P O Box 2541 Southport, QLD 4215
PH: 1300-786-454 F: 5591 8522
E: sales@quest surgical.com.au
ABN 33 118 090 618 ACN 118 090 618

CODE | DESCRIPTION | QTY.

SPEYES™ EYE SHIELDZ™

23-501	SPEyes™ Lenses	Box of 100
23-502	SPEyes™ Frames, Black	Box of 25
23-504	SPEyes™ Lens Dispenser Box Bracket	Each
23-505	Pre-assembled SPEyes™	Box of 100
23-506	Dispenser Box for pre-assembled SPEyes™	Box of 120
23-507	Dispenser Box Bracket for pre-assembled SPEyes™	Each



Specially designed, convenient dispenser box system makes lenses and frames easily accessible and disposable.

Product Number	UOM
DEF5526	Box of 10
FACE-SHIELD-FULL	Box of 300




WEARABLE CONSUMABLES

CLASSIC SURGICAL GOWNS (CONTINUED)

LEVEL 4 COOLMAX

Surgical Gown wrapped in a 91cm x 91cm SMS field with one huck towel. Sterile Level 4, trilaminate material.

High quality T-cut Surgical Gown ideal for high fluid procedures or lengthy procedures. Colour Coded neck band – blue for level 4 Gown.



Coolmax Surgical Gown Level 4

DEF5010A	Medium 110cm	25
DEF5011A	Large 130cm	25
DEF5012A	X-Large 140cm	25
DEF5013A	XL-XLong 160cm	25
DEF3198A	Coolmax Surgical Gown Pack – 1 Large and 1 X-Large Gown wrapped in a 200cm x 200cm SMS field with two huck towels	10

SPEyes EYEWEAR NON STERILE



Our popular DeRoyal SPEyes eyewear is designed to minimise exposure to pathogens from waste fluid. They combine innovation, quality, value, and comfort for clinicians where personal protective equipment is a vital part of self-care.

- Wrap-around disposable lens for complete protection
- UltraClear™ technology prevents glare, fog, and static
- Frame fits completely against forehead for splash protection from above
- Sealed lens mounting for added fluid protection
- Clip-On/Off design for quick and easy changing
- Coloured top strip allows finger-print-free handling
- Latex free



23-505

23-505	SPEyes Eye Shields pre-assembled assorted colours	100
--------	---	-----

SPEyes EYEWEAR NON STERILE



23-506

23-506	Dispenser box with pre-assembled SPEyes Eye Shields (30 per box 4 boxes per case)	120 total
23-501	SPEyes Eye Shield Replacement Lens	100
23-502	SPEyes Eye Shield Frame	25

FACESHIELD™



23-002

23-002	Face Mask Full with Elastic Strap	50
--------	-----------------------------------	----

GLASSES



DEF1755

DEF1755	Disposable Glasses Clear Polycarbonate Non Sterile	1
---------	--	---

BERET CAPS



DEF1142


DEF1142	Beret Cap 53cm (21") Blue Non Sterile	250
---------	---------------------------------------	-----

OVERSHOES NON STERILE



DEF1145

DEF1145	Non Skid Universal Overshoes Polypropylene	100
---------	--	-----



GOWNS

PATIENT GOWNS

Patient Gowns are lightweight, soft and comfortable and are available in a range of sizes to accommodate all body shapes. They are ideal for patients in hospitals, day procedures, GP clinics, radiology or any setting where a patient is required to wear a disposable gown.



PATIENT-Gown-SS-MED	Medium Short Sleeves	50
PATIENT-Gown-SS-LRG	Large Short Sleeves	50
PATIENT-Gown-SS-XL	X-Large Short Sleeves	50

SCRUB SUITS

Scrub Suits are made from a comfortable and soft fabric, and offer good protection. Ideal for nurses in the operating room, ER Departments, GP Clinics, Dentists, Vets and Students.



SCRUB-SUIT-SMLTOPS

SCRUB-SUIT-SML-TOPS	Top – Size Small	100
SCRUB-SUIT-MED-TOPS	Top – Size Medium	100
SCRUB-SUIT-LRG-TOPS	Top – Size Large	100
SCRUB-SUIT-XL-TOPS	Top – Size XL	100
SCRUB-SUIT-XXL-TOPS	Top – Size XXL	100



SCRUB-SUIT-SMLPANTS

SCRUB-SUIT-SML-PANTS	Pants – Size Small	100
SCRUB-SUIT-MED-PANTS	Pants – Size Medium	100
SCRUB-SUIT-LRG-PANTS	Pants – Size Large	100
SCRUB-SUIT-XL-PANTS	Pants – Size XL	100
SCRUB-SUIT-XXL-PANTS	Pants – Size XXL	100

SURGICAL GOWNS

Our single-use surgical gowns have specifically been developed for the Australian and New Zealand environment. Our range encompasses gowns for short and dry procedures through to high fluid and lengthy procedures.

OUR COMFORT RANGE – OUR SOFTEST AND MOST COMFORTABLE GOWNS

- Comfort Surgical Gown – Level 2.
- Comfort Unreinforced with Glove Assist – Level 3 gown with a reinforced sleeve to allow extra strength during gloving, but still cool, soft and ideal for short and dry procedures.
- Comfort Plus Fabric Reinforced Surgical Gown – Level 3 gown includes fluid resistant fabric reinforced panels yet remains lightweight and durable.
- Comfort Plus Poly Reinforced Surgical Gown – Level 3 gown which includes poly reinforced impervious panels for added protection.
- Comfort Ultimate Protection – Fully impervious high tech fabric, but still soft, comfortable and cool which meets AAMI Level 4 fabric standards.

OUR CLASSIC RANGE:

- Classic Standard Surgical Gown – suitable for most surgical procedures.
- Classic Plus Reinforced Surgical Gown – Includes a reinforced panel on the front and sleeves, for surgical cases requiring extra protection.
- Coolmax Surgical Gown – breathable and cool, these gowns are notable for their superb comfort and being fully impervious. With AAMI Level 4 fabric.

Designed in close co-operation with hospital staff, our Surgical Gowns are cool and breathable with an easy fit, while offering excellent protection against strikethrough from microbial and liquid penetration.



Beret Caps and Overshoes

- Single Use
- Non-Sterile and Latex Free



BERET CAPS

Product no.	Description	Quantity
DEF1142	53cm (21"), Blue	Box of 250



OVERSHOES

Product no.	Description	Quantity
DEF1145	Universal, Non Skid, Polypropylene	Box of 100