



CLINICAL DRAPES AND EQUIPMENT COVERS

BAGS

CLIP SEAL BAG



DEF640

DEF640	20cm x 25cm	100
DEF3286	20cm x 25cm, Packet of 20	50
DEF1703	28cm x 38cm	25

FREEZER BAG



DEF639

DEF639	30cm x 45cm	100
DEF66	30cm x 45cm with 2 x Rubber Bands #34	100
DEF1336	30cm x 45cm with 2 x Rubber Bands #34, Double Packed	100

DRAWSTRING BAG



DEF846

DEF846	25cm x 50cm	25
DEF864	44.5cm x 46.5cm	25
DEF861	44.5cm x 46.5cm	100
DEF1687	44.5cm x 46.5cm, Packet of 2	50
DEF848	50cm x 50cm	25
DEF847	50cm x 85cm	25

PLASTIC BAG



DEF1878

DEF1878	14cm x 30cm	100
DEF1276	14cm x 30cm with a Blank Adhesive Label	100
DEF638	28cm x 32.5cm	100
DEF650	Burns 30.5cm x 39cm	100
DEF642	45.5cm x 51cm with Cotton Tape 60cm Long	100
DEF808	Heavy Duty 57.5cm x 62cm	100

EQUIPMENT COVERS

BACK TABLE COVER

DEF5015	112cm x 228cm (44" x 90")	40
DEF3678	150cm x 300cm, Sapphire, HD	12
DEF5016AQ	152cm x 228cm (60" x 90"), Aqua	40
DEF5802	152cm x 228cm (60" x 90") Aqua Zoned Heavy Duty	20
DEF3544	150cm x 190cm, Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF3545	150cm x 230cm, Sapphire, Heavy Duty, Impervious and Absorbent	30
DEF5017	152cm x 228cm (60" x 90"), Heavy Duty	40
DEF5018	178cm x 280cm (70" x 110"), Heavy Duty	40

CAMERA & CORD COVER

		
DEF1164		
DEF1164	Camera Cord Cover 15cm x 230cm, Sealed End	50
DEF1201	Camera Cord Cover 15cm x 230cm with Cuff, Sealed End & Tapes, Double Packed	50

CHAIR COVER

		
DEF772		
DEF772	2 x 57.5cm x 62cm with Sealed Ends, with 2 x Rubber Bands #33	100

DRILL SLEEVE COVER

		
DEF887		
DEF1224	5cm x 80cm with 3 Labels	100
DEF1368	5cm x 110cm with Tape	100
DEF300	7.5cm x 75cm with Tape	100
DEF1516	7.5cm x 100cm with Adhesive both Ends	100
DEF1154	7.5cm x 150cm with 2 x Rubber Bands #31	100
DEF318	10cm x 60cm, Packet of 3	100
DEF887	10cm x 100cm	100
DEF889	10cm x 150cm	100
DEF618	10cm x 152cm with Tape, Double Packed	25


HANDLE SLEEVE COVER

		
DEF1053		
DEF1053	14cm x 30cm with Adhesive and Sealed End	100

MAYO STAND COVER

DEF5019AQ	61cm x 135cm with Sealed End	30
DEF5019	58.5cm x 142cm, Reinforced with Absorbent pad 50cm x 65cm, Cuff and Sealed End	50
DEF761	62cm x 91cm, Clear Plastic with Sealed End	100
DEF775	70cm x 137cm, Blue Plastic with Cuff and Sealed End	50
DEF3488	81cm x 165cm, Reinforced, Heavy Duty with Sealed End	25

PLASTIC SLEEVE COVER

		
DEF454		
DEF454	15cm x 100cm, with 2 x Rubber Bands #31, Double Packed	100
DEF1744	15cm x 200cm, with Rubber Band #31 and Sealed End	100

PROTECTIVE COVER

DEF613	15cm x 100cm, Plastic with Sealed End	50
DEF1744	15cm x 2m Sealed End, with 1 x Rubber Band #31, Double Packed	100
DEF763	104cm x 150cm, Protective	40
DEF1573	150cm x 208cm, Protective	50
DEF2582	27cm x 92cm, IV Pole Cover with Sealed End, Double Packed	100

TROLLEY COVER

		
DEF641		
DEF641	Trolley Bag Cover 61cm x 102cm with Sealed End	50
DEF1239	Shower Trolley Drape 2m x 3m	25



CUSTOMISED PROCEDURE DRAPES

THE FLEXIBILITY OF CHOICE

At Defries we have a comprehensive range of patient drapes which cover all procedures in the Operating Room and we can also design drapes in conjunction with Surgeons, Nurses and Theatre Managers. We provide our customers with superior products that match their clinical needs.

Our drapes are made to specific requirements such as size, fabric, fenestrations and extra details like velcro taping and clear plastic anaesthetic screens.

CHOICE OF FABRIC IS IMPORTANT

At Defries we use a variety of drape fabrics:

Azure drapes use a viscose impervious fabric with a polyethylene laminate, suitable for general surgery with minimal fluid or short procedures. This range includes the Azure minor procedure drapes.

Sapphire drapes use a bi-laminate or tri-laminate fabric with extra reinforcement. This range includes back table covers, mayo covers, some procedure drapes and bespoke designed drapes to meet specific customer requirements.

Aqua drapes are our newest drape range available and use an innovative drape fabric that is a tri-laminate – a three layer drape which offers an absorbent top level, an impervious middle coating and a soft and comfortable patient layer.

This high performance fabric is soft and drapeable with additional features such as reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. We have an excellent standard range of designs covering all procedures from general surgery and universal packs, to orthopaedics and specialised surgery such as neuro.

Supreme drapes are made from a tri-laminate fabric and can be tailored and customised to your surgeons' requirements. They are specifically designed for high level performance, coupled with great absorbency from fabric reinforcement around the fenestration, fluid containment pouches, positioning tape and incise area. The Supreme range includes drape sheets and specialty procedure drapes.

Our drape sheets, universal packs and specialty drapes include the instructions for correct aseptic application with markings and directions on the drape to assist ease of use.

All drape fabrics are conformable to patients' shapes. They are low linting and have excellent barrier properties to prevent strikethrough of airborne and blood borne bacteria.

Please contact your Defries representative for further information.

GENERAL DRAPES AND PACKS

AZURE MINOR PROCEDURE DRAPE

Soft viscose fabric combined with a polyethylene laminate to create a drape ideal for minor procedures.

DEF5094	38cm x 45cm, Azure Circumcision Drape with a 3cm Circular Fenestration	50
DEF5090	50cm x 60cm, Azure Minor Procedure Drape with a 6.5cm Circular Fenestration with Full Adhesive	25
DEF5091	91cm x 91cm, Azure Minor Procedure Drape with a 10cm x 15cm Oval Fenestration with Full Adhesive	40
DEF5092	200cm x 300cm, Azure Minor Procedure Drape with a 15cm x 25cm Fenestration with Full Adhesive	12

AQUA DRAPE SHEETS

Innovative tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.

DEF5024AQ	Half Drape, 100cm x 150cm	20
DEF5025AQ	Medium Drape, 100cm x 180cm	20
DEF5026AQ	3/4 Drape, 150cm x 200cm	20
DEF5028AQ	Large Drape, 183cm x 215cm	20
DEF5029AQ	XLarge Drape, 178cm x 254cm	12
DEF5030AQ	Top/Head Drape, 135cm x 250cm with Reinforcement and Adhesive	15
DEF5031AQ	Side Drape, 100cm x 130cm with Reinforcement and Adhesive	20
DEF5032AQ	Bottom/Foot Drape, 190cm x 185cm with Reinforcement and Adhesive	20

SUPREME DRAPE SHEETS

Supreme tri-laminate fabric specifically designed with excellent absorbency, low lint, drapeability and patient comfort.

DEF5024	Half Drape, 112cm x 145cm	20
DEF5025	Medium Drape, 112cm x 180cm	20
DEF5027	3/4 Drape, 142cm x 193cm with Reinforcement	20
DEF5028	Large Drape, 183cm x 215cm	20
DEF5030	Top/Head Drape, 135cm x 250cm with Reinforcement and Adhesive	30
DEF5031	Side Drape, 100cm x 110cm with Reinforcement and Adhesive	20

CLEAR PLASTIC DRAPE



DEF742

DEF742	45cm x 61cm	200
DEF741	61cm x 91cm	25
DEF079	61cm x 91cm	200
DEF749	70cm x 120cm	100
DEF745	76cm x 106cm	100
DEF3910	76cm x 152cm, Double Packed	100
DEF151	91cm x 115cm	100
DEF2430	Procedure Adhesive Drape, 91cm x 115cm with 2.5cm x 115cm adhesive down length of drape	50
DEF748	120cm x 140cm	100
DEF752	120cm x 140cm, Double Packed	100

UNIVERSAL DRAPE PACKS

UNIVERSAL PACK SHORT

DEF5212	<p>1 x Top Drape, 135cm x 250cm with Reinforcement 40cm x 70cm and Adhesive</p> <p>2 x Side Drapes, 100cm x 100cm with Reinforcement 40cm x 60cm and Adhesive</p> <p>1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement 40cm x 65cm and Adhesive</p> <p>1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement 50cm x 65cm</p> <p>Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement 80cm x 230cm</p>	10
---------	--	----

UNIVERSAL PACK LONG

DEF5213	<p>1 x Top Drape, 135cm x 250cm with Reinforcement 40cm x 70cm and Adhesive</p> <p>2 x Side Drapes, 100cm x 195cm with Reinforcement 40cm x 75cm and Adhesive</p> <p>1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement 40cm x 65cm and Adhesive</p> <p>1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement 50cm x 65cm</p> <p>Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement 75cm x 230cm</p>	5
---------	--	---

AQUA UNIVERSAL DRAPE PACK – LONG SIDES


DEF5400	<p>1 x Top Drape, 135cm x 250cm with Reinforcement 40cm x 70cm and Adhesive</p> <p>2 x Side Drapes, 100cm x 195cm with Reinforcement 40cm x 75cm and Adhesive</p> <p>1 x Bottom/Foot Drape, 185cm x 190cm with Reinforcement 40cm x 65cm and Adhesive</p> <p>1 x Mayo Stand Cover, 60cm x 128cm with Reinforcement 50cm x 65cm</p> <p>Wrapped in a Back Table Cover, 150cm x 230cm with Reinforcement 80cm x 230cm</p>	5
---------	--	---

EPIDURAL DRAPE

		
DEF1481		
DEF1481	76cm x 106cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100
DEF748FA-S	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	10
DEF3394	120cm x 140cm with 10cm x 15cm Fenestration and 2.5cm Adhesive around Fenestration and Corner Tabs	100

POLYBACKED DRAPE

An absorbent drape with a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor procedures.

		
DEF1146		
DEF1146	Polybacked Towel, 33cm x 46cm	25
DEF833	Polybacked Dressing Towel, 33cm x 46cm Large	750
DEF1287	Polybacked Towel, 40cm x 60cm	100
DEF1288	Polybacked Towel, 40cm x 60cm with a 10cm Circular Fenestration	100
DEF1289	Polybacked Towel, 40cm x 60cm with a 10cm x 15cm Fenestration and Adhesive on each corner of Drape 2.5cm x 2.5cm	100
DEF2072	Polybacked Procedural Drape, 75cm x 105cm with a 10cm x 15cm Fenestration and a .5cm Adhesive length along the 10cm sides of Fenestration	25

FABRIC SMS DRAPE

SMS fabric is a 3 layer polypropylene material bonded together to form a drape that is fluid repellent. Ideal for minor procedures.

		
DEF1594		
DEF1122	30.5cm x 61cm	100
DEF1594	61cm x 61cm	50

STOCKINETTE

		
DEF5060		
DEF5060	Medium 20cm x 122cm	25
DEF5061	Large 23cm x 122cm	25
DEF5062	XLarge 30cm x 122cm	25

LEGGINGS

DEF2481AQ	Aqua Lithotomy Leggings, 80cm x 140cm, Packet of 2	26
-----------	--	----

SUPPLEMENTARY

HAND DRYING AND STERILE FIELDS

BLUE HUCK TOWEL

High absorbency cotton towels.

		
DEF1300		
DEF1976	40cm x 68cm	30
DEF1300	40cm x 68cm	100
DEF2018	40cm x 68cm Blue, Packet of 5	40
DEF2219	40cm x 68cm with X-Ray, Blue	25

HAND DRYING AND STERILE FIELDS (CONTINUED)

GREEN HUCK TOWEL



DEF2619

DEF2619	40cm x 68cm, Packet of 2	100
DEF2501	40cm x 68cm, Packet of 4	50
DEF2140	40cm x 68cm, Packet of 5	40

WHITE HUCK TOWEL

DEF2753	40cm x 68cm	30
---------	-------------	----

DRESSING T TOWEL



DEF865

DEF865	32cm x 32.5cm	100
DEF243	32cm x 32.5cm	1000
DEF862	32cm x 32.5cm, Packet of 2	50
DEF879	32.5cm x 60cm Large	100
DEF857	32.5cm x 60cm Large, Packet of 2	50

DRESSING D TOWEL



DEF656

DEF656	32cm x 34cm, Packet of 2	50
DEF831	32cm x 34cm, Packet of 2	300
DEF832	32cm x 34cm, Packet of 10	100
DEF894	34cm x 60cm Large	50
DEF471	34cm x 60cm	500
DEF912	34cm x 60cm Large, Packet of 2	50
DEF910	34cm x 60cm Large, Packet of 5	25

DRESSING TOWEL

DEF470	40cm x 42cm, White	50
--------	--------------------	----



DEF470

we use different types of drape fabric

PE Film is a poly-ethylene material that is completely impervious and is used in our plastic drape range.

Poly backed Drapes are a 3 ply cotton top layer and a waterproof plastic backing. Soft and low linting, ideal for catheter and minor ward procedures.

SMS is an excellent all-purpose fabric that is cost-effective and ideal for simple procedures. It offers a sheer resistance with water repellent features.

Azure drapes use a viscose impervious fabric with a polyethylene laminate, suitable for short procedures with minimal fluid. This bi-laminate range includes the Azure minor procedure drapes.

Airlaid fabric features excellent high absorbency and is designed specifically for the Angiography and Cath Lab procedures.

Sapphire fabric is chosen for our back table covers, mayo stand covers and some specialized drapes, this is a versatile super-absorbent fabric which provides excellent barrier protection.

Aqua drapes are our newest drape range available and use a bi-laminate or tri-laminate fabric with extra reinforcement. The tri-laminate fabric is soft to touch, whilst offering the superior protection of great absorbency from fabric reinforcement, fluid containment built-in pouches, positioning tape and incise area. The bi-laminate range includes back table covers & mayo covers, whilst the tri-laminate range are our drape sheets, universal packs and specialty procedure drapes. We can also offer bespoke designed drapes to meet specific customer requirements using this fabric.

Aqua HD is a superior version of the Aqua - designed for specialized drapes and Universal Packs, and offers full protection whilst still cost-effective.

Supreme drapes are our ultimate tri-laminate fabric. They are specifically designed for high-level performance and offer a high level of protection. The Supreme range includes drape sheets and specialty procedure drapes.