

Quest Surgical

Surgical Hand Immobilizers

with anti-tilt Tab Lead-free and Latex-free
Quest Surgical has TGA Registration
and the BEST pricing for
the Alumi and Dura Hands.



Alumi-Hand is sterile and ready to use

With the affordable Alumi-Hand you trade the hazards of lead for a sterile, disposable, X-ray compatible hand immobilizer. The Alumi-Hand is made with a pliable aluminium interior for easy hand positioning and a foam exterior for patient comfort. Available in large



TGA for Alumi Hand: ARTG198376

Dura-Hand – Easy to sterilize and is reusable

Just like the Alumi-Hand, the Dura-Hand is lead-free and X-ray compatible. The Dura-Hand is made of pliable aluminium for easy positioning. Available in Pediatric, Large and Extra-Large sizes.

TGA for Dura Hand: ARTG199569

Product Code	Description	Qty per Box
QS-760	Large Alumi-Hands 9 1/2 Inch x 8 1/2 inch Sterile Single Use	4
DH-100	Pedi DuraHands 7 1/4 Inch x 6 3/4 inch Reusable	1
DH-225	Large DuraHands 9 1/2 Inch x 8 1/2 inch Reusable	1
DH-350	Extra-Large DuraHands 13 Inch x 10 1/2 inch Reusable	1

P O Box 2541 Southport QLD 4215
 Ph: 1300 786 454 Fax: 5591 8522
 Email: sales@questsurgical.com.au

ABN 33 118 090 618 ACN 118 090 618



Quality
 ISO 9001
 SAI GLOBAL